

Product Datasheet

Rabbit APOA1 ELISA Kit (orb437993)

| | |
|--------------------------|------------------------------------------------------------------------------------------|
| Description | Rabbit Apolipoprotein A1 (APOA1) competitive inhibition ELISA Kit |
| Reactivity | Rabbit |
| Range | 31.25-2000ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | P09809 |
| Sensitivity | 9.35ng/mL |
| Entrez | 100009253 |
| Expiration Date | Please enquire. |