

Product Datasheet

Rat APOC2 ELISA Kit (orb438020)

Description	Rat Apolipoprotein C2 (APOC2) Sandwich ELISA Kit
Reactivity	Rat
Range	39.062-2500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	G3V8D4
Sensitivity	14pg/mL
Entrez	292697
Expiration Date	Please enquire.