

Product Datasheet

Rat AQP2 ELISA Kit (orb438051)

Description	Rat Aquaporin 2, Collecting Duct (AQP2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P34080
Sensitivity	0.235ng/mL
Entrez	25386
Expiration Date	Please enquire.