

Product Datasheet

Rat bEP ELISA Kit (orb438185)

Description	Rat Beta-Endorphin (bEP) competitive inhibition ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernatants and other biological fluids.
Sensitivity	4.35pg/mL
Expiration Date	Please enquire.