

Product Datasheet

Mouse Bid ELISA Kit (orb438191)

Description	Mouse BH3 Interacting Domain Death Agonist (Bid) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P70444
Sensitivity	0.112ng/mL
Entrez	12122
Expiration Date	Please enquire.