

## Product Datasheet

### Porcine C5a ELISA Kit (orb438295)

<b>Description</b>	Porcine Complement Component 5a (C5a) Sandwich ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3.12-200ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Sensitivity</b>	0.072ng/mL
<b>Expiration Date</b>	Please enquire.