

## Product Datasheet

### Human CA19-5 ELISA Kit (orb438315)

<b>Description</b>	Human Carbohydrate Antigen 19-5 (CA19-5) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.125-200ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Sensitivity</b>	0.86ng/mL
<b>Expiration Date</b>	Please enquire.