

Product Datasheet

Rat CASP8 ELISA Kit (orb438398)

Description	Rat Caspase 8 (CASP8) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9JHX4
Sensitivity	0.053ng/mL
Entrez	64044
Expiration Date	Please enquire.