

Product Datasheet

Human CD200R2 ELISA Kit (orb438488)

Description	Human CD200 Receptor 2 (CD200R2) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q6Q8B3
Sensitivity	0.34ng/mL
Entrez	344807
Expiration Date	Please enquire.