

## Product Datasheet

### Human CD5L ELISA Kit (orb438509)

<b>Description</b>	Human CD5 Antigen Like Protein (CD5L) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>O43866</b>
<b>Sensitivity</b>	28.6pg/mL
<b>Entrez</b>	<b>922</b>
<b>Expiration Date</b>	Please enquire.