

Product Datasheet

Human CEACAM8 ELISA Kit (orb438560)

Description	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 8 (CEACAM8) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P31997
Sensitivity	5.9pg/mL
Entrez	1088
Expiration Date	Please enquire.