

Product Datasheet

Rat CES1 ELISA Kit (orb438575)

| | |
|--------------------------|--|
| Description | Rat Carboxylesterase 1 (CES1) Sandwich ELISA Kit |
| Reactivity | Rat |
| Range | 0.312-20ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | Q63108 |
| Sensitivity | 0.107ng/mL |
| Entrez | 29225 |
| Expiration Date | Please enquire. |