

## Product Datasheet

### Mouse CHEM ELISA Kit (orb438616)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Mouse Chemerin (CHEM) Sandwich ELISA Kit  |
| <b>Reactivity</b>        | Mouse   |
| <b>Range</b>             | 78.125-5000pg/mL  |
| <b>Form/Appearance</b>   | Traditional   |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.                |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| <b>Uniprot ID</b>        | <a href="#">Q9DD06</a>  |
| <b>Sensitivity</b>       | 26pg/mL   |
| <b>Entrez</b>            | <a href="#">71660</a>   |
| <b>Expiration Date</b>   | Please enquire.   |