

Product Datasheet

Rat ChT ELISA Kit (orb438641)

Description	Rat High Affinity Choline Transporter (ChT) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9GZV3
Sensitivity	28pg/mL
Entrez	60482
Expiration Date	Please enquire.