

Product Datasheet

Mouse CLDN3 ELISA Kit (orb438671)

Description	Mouse Claudin 3 (CLDN3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9Z0G9
Sensitivity	0.061ng/mL
Entrez	12739
Expiration Date	Please enquire.