

Product Datasheet

Canine CLU ELISA Kit (orb438690)

Description	Canine Clusterin (CLU) Sandwich ELISA Kit
Reactivity	Canine
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P25473
Sensitivity	0.5ng/mL
Entrez	442971
Expiration Date	Please enquire.