

Product Datasheet

Mouse CLU ELISA Kit (orb438692)

Description	Mouse Clusterin (CLU) Sandwich ELISA Kit
Reactivity	Mouse
Range	18.75-1200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine, breast milk, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, urine, breast milk, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernatants and other biological fluids.
Uniprot ID	Q06890
Sensitivity	6.1ng/mL
Entrez	12759
Expiration Date	Please enquire.