

Product Datasheet

Human CNDP1 ELISA Kit (orb438697)

Description	Human Carnosine Dipeptidase 1 (CNDP1) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cerebrospinal fluid and other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid and other biological fluids.
Uniprot ID	Q96KN2
Sensitivity	1.45ng/mL
Entrez	84735
Expiration Date	Please enquire.