


Product Datasheet

HRP antibody (orb43872)

Description	Mouse Monoclonal to HRP.
Reactivity	All, Plant
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Horseradish peroxidase (Grade I)
Target	HRP
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunocytochemistry: The antibody HP-03 was successfully used for staining of formaldehyde-fixed, Triton-permeabilized cells transfected with HRP gene.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Dilution Range	Immunocytochemistry: The antibody HP-03 was successfully used for staining of formaldehyde-fixed, Triton-permeabilized cells transfected with HRP gene.
Expiration Date	12 months from date of receipt.



Western blotting analysis of purified ho...