

Product Datasheet

Human COL15 ELISA Kit (orb438726)

| | |
|--------------------------|--|
| Description | Human Collagen Type XV (COL15) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | P39059 |
| Sensitivity | 0.121ng/mL |
| Entrez | 1306 |
| Expiration Date | Please enquire. |