

Product Datasheet

Human COQ10B ELISA Kit (orb438783)

Description	Human Coenzyme Q10 Homolog B (COQ10B) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9H8M1
Sensitivity	0.124ng/mL
Entrez	80219
Expiration Date	Please enquire.