

Product Datasheet

Human CPLX2 ELISA Kit (orb438806)

Description	Human Complexin 2 (CPLX2) competitive inhibition ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6PUV4
Sensitivity	0.084ng/mL
Entrez	10814
Expiration Date	Please enquire.