

## Product Datasheet

### Human CRHR1 ELISA Kit (orb438843)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human Corticotropin Releasing Hormone Receptor 1 (CRHR1) Sandwich ELISA Kit              |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 0.156-10ng/mL  |
| <b>Form/Appearance</b>   | Traditional  |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | tissue homogenates, cell lysates and other biological fluids.                            |
| <b>Sample Types</b>      | tissue homogenates, cell lysates and other biological fluids.                            |
| <b>Uniprot ID</b>        | <a href="#">P34998</a>   |
| <b>Sensitivity</b>       | 0.060ng/mL   |
| <b>Entrez</b>            | <a href="#">1394</a>   |
| <b>Expiration Date</b>   | Please enquire.  |