

Product Datasheet

Human CRT ELISA Kit (orb438873)

Description	Human Calreticulin (CRT) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P27797
Sensitivity	0.061ng/mL
Entrez	811
Expiration Date	Please enquire.