

Product Datasheet

Porcine CST3 ELISA Kit (orb438895)

Description	Porcine Cystatin 3 (CST3) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	Q0Z8R0
Sensitivity	0.056ng/mL
Entrez	733672
Expiration Date	Please enquire.