

Product Datasheet

Human CSTB ELISA Kit (orb438901)

Description	Human Cystatin B (CSTB) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P04080
Sensitivity	0.118ng/mL
Entrez	1476
Expiration Date	Please enquire.