

Product Datasheet

Human CTXIII ELISA Kit (orb438951)

Description	Human Cross Linked C-Telopeptide Of Type III Collagen (CTXIII) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.043ng/mL
Expiration Date	Please enquire.