

Product Datasheet

Mouse CTXII ELISA Kit (orb438952)

Description	Mouse Cross Linked C-Telopeptide Of Type II Collagen (CTXII) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	4.26pg/mL
Expiration Date	Please enquire.