

Product Datasheet

Rat CTXII ELISA Kit (orb438953)

Description	Rat Cross Linked C-Telopeptide Of Type II Collagen (CTXII) Sandwich ELISA Kit
Reactivity	Rat
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	33.5pg/mL
Expiration Date	Please enquire.