

Product Datasheet

Human CYP1A1 ELISA Kit (orb438997)

Description	Human Cytochrome P450 1A1 (CYP1A1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P04798
Sensitivity	0.056ng/mL
Entrez	1543
Expiration Date	Please enquire.