

## Product Datasheet

### Rat CYP1A1 ELISA Kit (orb438998)

<b>Description</b>	Rat Cytochrome P450 1A1 (CYP1A1) Sandwich ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P00185</a>
<b>Sensitivity</b>	0.129ng/mL
<b>Entrez</b>	<a href="#">24296</a>
<b>Expiration Date</b>	Please enquire.