

Product Datasheet

Rat DBP ELISA Kit (orb439059)

Description	Rat D Site Of Albumin Promoter Binding Protein (DBP) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P16443
Sensitivity	0.105ng/mL
Entrez	24309
Expiration Date	Please enquire.