

Product Datasheet

Human DEFa6 ELISA Kit (orb439087)

Description	Human Defensin Alpha 6, Paneth Cell Specific (DEFa6) Sandwich ELISA Kit
Reactivity	Human
Range	2.343-150ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q01524
Sensitivity	1.07ng/mL
Entrez	1671
Expiration Date	Please enquire.