

Product Datasheet

Human DEFb1 ELISA Kit (orb439096)

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| Description | Human Defensin Beta 1 (DEFb1) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 62.5-4000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | P60022 |
| Sensitivity | 25.3pg/mL |
| Entrez | 1672 |
| Expiration Date | Please enquire. |