

Product Datasheet

Human DPAGT1 ELISA Kit (orb439176)

| | |
|--------------------------|--|
| Description | Human Dolichyl Phosphate-N-Acetylglucosaminephosphotransferase 1 (DPAGT1) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | tissue homogenates, cell lysates and other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | Q9H3H5 |
| Sensitivity | 0.059ng/mL |
| Entrez | 1798 |
| Expiration Date | Please enquire. |