

Product Datasheet

Human ECE1 ELISA Kit (orb439240)

Description	Human Endothelin Converting Enzyme 1 (ECE1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P42892
Sensitivity	11.3pg/mL
Entrez	1889
Expiration Date	Please enquire.