

Product Datasheet

Mouse EGFL7 ELISA Kit (orb439263)

Description	Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9QXT5
Sensitivity	15.5pg/mL
Entrez	353156
Expiration Date	Please enquire.