

Product Datasheet

Canine ELA2 ELISA Kit (orb439275)

Description	Canine Elastase 2, Neutrophil (ELA2) Sandwich ELISA Kit
Reactivity	Canine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8MJD1
Sensitivity	1.23ng/mL
Entrez	442980
Expiration Date	Please enquire.