

Product Datasheet

Human EPHA10 ELISA Kit (orb439322)

Description	Human Ephrin Type A Receptor 10 (EPHA10) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q5JZY3
Sensitivity	0.109ng/mL
Entrez	284656
Expiration Date	Please enquire.