

Product Datasheet

Canine EPO ELISA Kit (orb439329)

Description	Canine Erythropoietin (EPO) Sandwich ELISA Kit
Reactivity	Canine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P33707
Sensitivity	12.8pg/mL
Entrez	404002
Expiration Date	Please enquire.