


Product Datasheet

CD20 antibody (Biotin) (orb43938)

Description	Mouse Monoclonal to CD20.
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Normal human lymphocytes from lymph node.
Target	CD20
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P11836
Entrez	931
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface stainin...