

Product Datasheet

Mouse FABP6 ELISA Kit (orb439434)

Description	Mouse Fatty Acid Binding Protein 6, Ileal (FABP6) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P51162
Sensitivity	23pg/mL
Entrez	16204
Expiration Date	Please enquire.