

Product Datasheet

Mouse FBLN3 ELISA Kit (orb439476)

Description	Mouse Fibulin 3 (FBLN3) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	Q8BPB5
Sensitivity	6.6pg/mL
Entrez	216616
Expiration Date	Please enquire.