

Product Datasheet

Human FKBP1A ELISA Kit (orb439533)

Description	Human FK506 Binding Protein 1A (FKBP1A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P62942
Sensitivity	0.109ng/mL
Entrez	2280
Expiration Date	Please enquire.