

Product Datasheet

Human FLII ELISA Kit (orb439541)

Description	Human Flightless I Homolog (FLII) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q13045
Sensitivity	29.2pg/mL
Entrez	2314
Expiration Date	Please enquire.