

Product Datasheet

Mouse FOLR1 ELISA Kit (orb439564)

Description	Mouse Folate Receptor 1, Adult (FOLR1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35846
Sensitivity	31pg/mL
Entrez	14275
Expiration Date	Please enquire.