

Product Datasheet

Mouse FSH ELISA Kit (orb439592)

| | |
|--------------------------|--|
| Description | Mouse Follicle Stimulating Hormone (FSH) competitive inhibition ELISA Kit |
| Reactivity | Mouse |
| Range | 3.125-200ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | Q60687 |
| Sensitivity | 1.49ng/mL |
| Entrez | 14308 |
| Expiration Date | Please enquire. |