

Product Datasheet

Human FSHR ELISA Kit (orb439595)

Description	Human Follicle Stimulating Hormone Receptor (FSHR) competitive inhibition ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P23945
Sensitivity	0.124ng/mL
Entrez	2482
Expiration Date	Please enquire.