

Product Datasheet

Human FTH ELISA Kit (orb439604)

Description	Human Ferritin, Heavy Polypeptide (FTH) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02794
Sensitivity	12.9pg/mL
Entrez	2495
Expiration Date	Please enquire.