

Product Datasheet

Human gGT5 ELISA Kit (orb439747)

Description	Human Gamma-Glutamyltransferase 5 (gGT5) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P36269
Sensitivity	0.55ng/mL
Entrez	2687
Expiration Date	Please enquire.